

DATA LOGGER



EliTech: How to Configure

ElitechLog V4.0.0

Download Summary **Parameter** Graph Database Import Email System About us Help Exit

Logger Clock: 2019-08-23 16:44:07

Log Interval: 0 H 5 M 0 S

Press Button: Enable

Temperature Unit: °C

Tone setting: Alarm tone

Alarm Threshold: High Limit(Temperature): 8.0

S/N: EFE18AD08122

Run Time: 111D 2H 35M 0S

Start Delay(H): 0.5

Storage Mode

Audible alarm interval

Low Limit(Temperature): 2.0

STEPS	DIRECTIONS
Step 1	Click: Parameter
Step 2	Change Log Interval to 5 minutes (0H 5M 0S)
Step 3	Press Button: Change selection to Enable
Step 4	Start Delay (H): Choose either 0.0 or 0.5 (equivalent to 30minutes)*
Step 5	Temperature Unit: Change selection to Celsius (°C)
Step 6	Click Configure to save
Step 7	High Limit (Temperature): change to 8°C
Step 8	Low Limit (Temperature): change to 2°C
	<i>These configurations will be saved for future – only need set parameters ONCE</i>

* A delay of 30 minutes allows for the logger to cool down once placed inside fridge.